
OCULOPLASTIK

(en)	Symblypharon Ring.....	1
(fr)	Anneau de symblépharon	5
(de)	Symblypharon-Ring	10
(es)	Anillo de simbléfaron.....	14
(it)	Anello per simblefaron	18
(nl)	Symblypharon Ring.....	23
(pt)	Anel para simbléfaro	27

(en) SYMBLEPHARON RING

INSTRUCTIONS FOR USE

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

CAUTION, consult accompanying documents.

This Product is sold Non Sterile.

oculoplastik.com

Get our instructions for use at oculoplastik.com.

Do not autoclave.

Do not reuse.

R_xonly For prescription use only.

DESCRIPTION

The Symblepharon Ring is a medical device made of acrylic (PMMA). The instruments' family has various sizes (small, medium, large).

APPEARANCE

Symblepharon Ring is sold non sterile. It should not be cleaned prior to sterilization. The shelf life is permanent.

INDICATIONS

Symblepharon Ring is intended to prevent retraction of the conjunctival cul-de-sacs without touching the cornea.

MODE OF ACTION

The ring serves as a barrier between the eyelid and the ocular globe allowing the patient to see.

CONTRAINDICATIONS

Cold sterilization is not recommended, because chemical ocular burns have been reported, due to inadequate rinsing or leaching out of sterilizing solution. It also may damage plastic.

Do not use alcohol, as it will damage acrylic (PMMA).

Do not use any device that is scratched or damaged.

ADVERSE EFFECTS

Potential complications associated to this device include, but are not limited to the following:

- A corneal abrasion may be caused by a scratched device.

Unexpected side effects and complications related to this device must be reported to Oculo-Plastik, Inc.

RECOMMENDATIONS FOR USE

STERILIZATION

The Non Sterile device should remain in the original sealed packaging until it will be STERILIZED. It should not be cleaned prior to sterilization. It should be sterilized alone in a separate Peel Pouch that does not contain other devices.

- The device can only be sterilized one time.
- Any device that has become contaminated by blood, tissue and/or bodily fluids should be thrown out and should not be processed in any way.
- The device should not be processed or transported or come in contact with any type of contaminated devices and/or equipment.
- The device should be transported separate from soiled devices in its original sealed packaging in order to avoid contamination.
- The device should remain in the original packaging, which is sealed until it will be STERILIZED in order to avoid contamination.
- The sealed package should be opened under clean conditions.
- Care must be taken to prevent possible contamination by talc, dust, lint, surgical glove powder, drape and sponge lint, fingerprints, skin oils, and other surface contaminants that may be deposited during handling prior to sterilization.
- Talc, dust, lint, surgical glove powder, drape and sponge lint, fingerprints, skin oils and other surface contaminants deposited by improper handling may adversely affect the suitability of the device.
- Strict adherence to clean, aseptic techniques should be maintained to prevent contamination of the device and possible complications.

- ONCE THE DEVICE HAS BEEN REMOVED FROM THE ORIGINAL PACKAGING IT SHOULD BE STERILIZED IMMEDIATELY in a SEPARATE peel pouch that does not contain other devices.

ETHYLENE OXIDE

To sterilize pouched devices, place them on the sterilizer rack, in a single layer, as pressure from other devices may damage them.

Article configuration: Individually pouched.

In Chamber Conditioning:

Temperature:	55°C
Relative Humidity:	80%
Vacuum Set Point:	1.3 psia (9 kPa)
Preconditioning Time:	30 minutes

Exposure:

Temperature:	55°C
Relative Humidity:	80%
Ethylene Oxide Concentration:	740 mg/L
Gas Exposure Time:	60 minutes
Gas Type:	100% EO

Aeration:

Aeration Time:	24 hours
Aeration Temperature:	55 ± 4°C

STERRAD

Our recommendation with Sterrad is a ONE TIME sterilization. (If sterilized more than once the biomaterial will become damaged and should not be used.) The device must be sterilized with a full cycle of the STERRAD 100S. The full cycle of the STERRAD 100S Sterilizer comprises TWO identical exposure phases (6 minute injection, 2 minute diffusion, and 2 minute plasma) preceded by a 10 minute pre-exposure plasma phase (1440 µL of 59% hydrogen peroxide). Follow the manufacturer's sterilization recommendations with the Sterrad 100S Sterilization System.

Insertion technique

The physician must insert the Symblepharon Ring between the eye and the eyelid. The ring must not touch the cornea.

PRECAUTIONS FOR USE

Do not use any device that has been dropped accidentally without thorough examination of its condition.

Symblepharon ring should not be left in place more than 30 days continuously.

The Symblepharon Ring is MR Safe and does not pose any known hazard in a magnetic resonance environment.

Symblepharon Ring is a **single use** medical device for a **single patient**:

- A second use of Symblepharon Ring at a later time may cause an infection due to microbial contamination.

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ATTENTION, consulter les documents joints.

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