
OCULOPLASTIK


(en)	Acrylic Sphere	1
(fr)	Sphère en acrylique.....	6
(de)	Acryl-Kugelimplantat.....	11
(es)	Esfera acrílica.....	16
(it)	Sfera acrilica.....	21
(nl)	Acryl bol	26
(pt)	Esfera acrílica.....	31

(en) ACRYLIC SPHERE

INSTRUCTIONS FOR USE

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

CAUTION, consult accompanying documents.

This Product is sold Non Sterile.

oculoplastik.com

Get our instructions for use at oculoplastik.com.

Do not autoclave.

Do not reuse.

Rx only For prescription use only.

DESCRIPTION

The acrylic sphere (PMMA) is a non-integrated ocular implant used for enucleation, evisceration or as a secondary implant. It exists in a range of different sizes: 12mm (1.1cc), 14mm (1.5cc), 16mm (2.1cc), 18mm (2.8cc), 20mm (3.6cc) and 22mm (5.7cc).

APPEARANCE

The sphere is delivered non sterile. It should not be cleaned prior to sterilization.

INDICATIONS

The sphere is intended to be used for enucleation, evisceration or as a secondary implant. It aims to replace the lost volume during surgery and to give motility to the ocular prosthesis. In areas of the world where ocular prosthesis are not usually custom-made, it is recommended to use a sphere.

MODE OF ACTION

The sphere is the most commonly used ocular implant for this type of surgery. But the motility of

the prosthesis is more or less poor. Some implantation techniques allow the motility to the prosthesis.

CONTRAINDICATIONS

Some surgeons imbricate ocular muscles over the sphere hoping to give motility to the sphere, but the sphere may migrate. A ptosis to the upper eyelid with a migrated sphere (from the imbrication of the ocular muscles on the sphere) constitute a real syndrome. To place the muscles over the PMMA implant, the surgeon should suture the muscles to a wrap (such as mesh, sclera, fascia lata, dura, etc.) covering the sphere.

ADVERSE EFFECTS

If the implantation goes well, there will generally be few adverse effects. It is suggested that the patient visits his ophthalmologist and his ocularist yearly to minimize possible complications. Potential adverse effects for an implantation of this implant include but are not limited to:

- Tissue necrosis
- Tissue dehiscence or erosion
- Implant exposure

- Implant extrusion
- Deep superior sulcus deformity
- Ptosis of the superior and inferior eyelid
- Poor motility
- Migration or displacement of the implant
- Volume deficiencies
- Socket cul-de-sacs contraction
- Discomfort
- Infection

The smooth spheres have tendency to migrate or to displace with time. Implant's migration or displacement comes when a muscle disinserts (detaches itself) or when a sphere is not covered by an envelope.

Unexpected side effects and complications related to this medical device must be reported to Oculo-Plastik, Inc.

RECOMMENDATIONS FOR USE

STERILIZATION

The Non Sterile device should remain in the original sealed packaging until it will be STERILIZED. It should not be cleaned prior to sterilization. It

should be sterilized alone in a separate Peel Pouch that does not contain other devices.

- The device can only be sterilized one time.
- Any device that has become contaminated by blood, tissue and/or bodily fluids should be thrown out and should not be processed in any way.
- The device should not be processed or transported or come in contact with any type of contaminated devices and/or equipment.
- The device should be transported separate from soiled devices in its original sealed packaging in order to avoid contamination.
- The device should remain in the original packaging, which is sealed until it will be STERILIZED in order to avoid contamination.
- The sealed package should be opened under clean conditions.
- Care must be taken to prevent possible contamination by talc, dust, lint, surgical glove powder, drape and sponge lint, fingerprints, skin oils, and other surface contaminants that may be deposited during handling prior to sterilization.
- Talc, dust, lint, surgical glove powder, drape and sponge lint, fingerprints, skin oils and

other surface contaminants deposited by improper handling may adversely affect the suitability of the device.

- Strict adherence to clean, aseptic techniques should be maintained to prevent contamination of the device and possible complications.
- ONCE THE DEVICE HAS BEEN REMOVED FROM THE ORIGINAL PACKAGING IT SHOULD BE STERILIZED IMMEDIATELY in a SEPARATE peel pouch that does not contain other devices.

ETHYLENE OXIDE

To sterilize pouched devices, place them on the sterilizer rack, in a single layer, as pressure from other devices may damage them.

Article configuration: Individually pouched.

In Chamber Conditioning:

Temperature: 55°C
Relative Humidity: 80%
Vacuum Set Point: 1.3 psia (9 kPa)
Preconditioning Time: 30 minutes

Exposure:

Temperature: 55°C
Relative Humidity: 80%
Ethylene Oxide Concentration: 740 mg/L
Gas Exposure Time: 60 minutes
Gas Type: 100% EO

Aeration:

Aeration Time: 24 hours
Aeration Temperature: 55 ± 4°C

STERRAD

Our recommendation with Sterrad is a ONE TIME sterilization. (If sterilized more than once the bio-material will become damaged and should not be used.) The device must be sterilized with a full cycle of the STERRAD 100S. The full cycle of the STERRAD 100S Sterilizer comprises TWO identical exposure phases (6 minute injection, 2 minute diffusion, and 2 minute plasma) preceded by a 10 minute pre-exposure plasma phase (1440 µL of 59% hydrogen peroxide). Follow the manufacturer's sterilization recommendations with the Sterrad 100S Sterilization System.

Implantation technique

The enucleation or evisceration technique is standard.

The appropriate size of the sphere must be chosen by the surgeon. Because of the orbital tissue contraction, and in order to avoid a dehiscence of closure, the size of the sphere must be adapted to the patient anatomy. Too small an implant may not occupy enough space of the orbital socket; too large an implant may extrude. The choice should allow a closure without tension. If the sphere is wrapped, the thickness of the wrap must be added to the sphere diameter in order to calculate the wanted volume. In all cases, the size must be in proportion with the volume to occupy.

The sphere can be implanted non-wrapped, especially for evisceration.

For evisceration, it is generally not recommended to preserve cornea, especially when there is a history of corneal ulcers.

For enucleation, you can wrap the sphere with absorbable mesh or autogenous material (sclera, fascia-lata, temporal fascia, dura).

Use an acrylic bilateral conformer 22 or 18/20 after the surgery. The use of a silicone conformer is never recommended.

PRECAUTIONS FOR USE

This implant is sold non sterile. Ethylene oxide sterilization (EO) is recommended.

The Acrylic Sphere is MR Safe and does not pose any known hazard in a magnetic resonance environment.

The Acrylic Sphere is a **single use** medical device for a **single patient**:

- A second use of Acrylic Sphere at a later time may cause an infection due to microbial contamination.

The storage's duration is permanent.

Do not use any device that has been dropped accidentally without thorough examination of its condition. Do not use any device that is scratched or damaged.

Cold sterilization is not recommended for a bio-material intended for implantation.

The patient should visit yearly his ophthalmologist and his ocularist in order to minimize subsequent complications. If there are any changes in the performance of the implant, the patient must communicate with one or another of those specialists.

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ATTENTION, consulter les documents joints.

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