
OCULOPLASTIK

(en)	Codère-Durette Orbital Floor Implant.....	1
(fr)	Implant de plancher d'orbite Codère-Durette.	6
(de)	Codère-Durette-Orbitabodenimplantat.....	11
(es)	Implante del piso orbital Codère-Durette.....	16
(it)	Impianto del pavimento orbitario Codère-Durette.....	21
(nl)	Oogkasbodemimplantaat Van Codère-Durette.....	26
(pt)	Implante de pavimento orbitário Codère-Durette.....	31

(en) CODÈRE-DURETTE ORBITAL FLOOR IMPLANT

INSTRUCTIONS FOR USE

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

CAUTION, consult accompanying documents.

This Product is sold Non Sterile.

oculoplastik.com

Get our instructions for use at oculoplastik.com.

Do not autoclave.

Do not reuse.

For prescription use only.

DESCRIPTION

These implants are prefabricated in methylmethacrylate (PMMA). They are unilateral; placement is either left of right, not both. They were designed with the orbital floors of an adult skull reproduction.

They are easily inserted as they fit snugly on the orbital floor and, when well fixated, minimize movement and possible protrusion.

They are available with three (3) volumes (for each side). They come in six (6) models; small - left, small – right, medium – left, medium – right, large - left, large – right. The added volume is situated on the posterior-superior aspects of the implants. The face under and the contours are identical for all volumes. They have a smooth surface and are permanent.

It is suggested to fix it with a titanium screw, as used in cranio-facial surgery of the bones. The holes are 3/32" (approximately 2,4 mm). Use an appropriate size screw and one of 7 mm in length.

APPEARANCE

This implant is sold non sterile. It should not be cleaned prior to sterilization.

INDICATIONS

Used only in cases of enucleated, eviscerated or non-seeing eyes. These implants are for subperiosteal placement. They serve to fill the superior sulcus depressions. They are available with three (3) volumes (for each side) to provide anterior displacement as well as elevation of the orbital contents. The implants are permanent, although a surgeon may elect to remove one or replace it with another, as it is a smooth implant.

MODE OF ACTION

The Codère-Durette® should improve symmetry of the eyes of the patient (cosmetically) by filling the volume deficit when a superior sulcus deformity is present.

CONTRAINDICATIONS

The Codère-Durette orbital floor implant should not be used for a seeing eye (like in the case of an orbital floor fracture), as it may cause diplopia. Also, a patient with shallow inferior cul de sacs may worsen his cul de sacs condition after implantation

ADVERSE EFFECTS

If the implantation goes well, there will generally be no adverse effect. It is suggested that the patient (like any patient with an artificial eye) visits his ophthalmologist and his ocularist yearly to prevent possible complications. Potential adverse effects for an implantation of this implant include but are not limited to:

- Loss of sensation of the cheek possibly because the implant touches a nerve at floor of the orbit
- Small cicatrix within the ridule of the lower lid
- Volume deficiencies due to over or under correction
- Migration or displacement of the implant
- Infection
- Difficulty in fitting the custom ocular prosthesis

Possible causes of these complications are: inappropriate implantation or fixation technique, insufficient inferior cul de sacs, too small an implant was chosen.

Unexpected side effects and complications related to this medical device must be reported to Oculo-Plastik, Inc.

RECOMMENDATIONS FOR USE

STERILIZATION

The Non Sterile device should remain in the original sealed packaging until it will be STERILIZED. It should not be cleaned prior to sterilization. It should be sterilized alone in a separate Peel Pouch that does not contain other devices.

- The device can only be sterilized one time.
- Any device that has become contaminated by blood, tissue and/or bodily fluids should be thrown out and should not be processed in any way.
- The device should not be processed or transported or come in contact with any type of contaminated devices and/or equipment.
- The device should be transported separate from soiled devices in its original sealed packaging in order to avoid contamination.
- The device should remain in the original packaging, which is sealed until it will be STERILIZED in order to avoid contamination.

- The sealed package should be opened under clean conditions.
- Care must be taken to prevent possible contamination by talc, dust, lint, surgical glove powder, drape and sponge lint, fingerprints, skin oils, and other surface contaminants that may be deposited during handling prior to sterilization.
- Talc, dust, lint, surgical glove powder, drape and sponge lint, fingerprints, skin oils and other surface contaminants deposited by improper handling may adversely affect the suitability of the device.
- Strict adherence to clean, aseptic techniques should be maintained to prevent contamination of the device and possible complications.
- ONCE THE DEVICE HAS BEEN REMOVED FROM THE ORIGINAL PACKAGING IT SHOULD BE STERILIZED IMMEDIATELY in a SEPARATE peel pouch that does not contain other devices.

ETHYLENE OXIDE

To sterilize pouched devices, place them on the sterilizer rack, in a single layer, as pressure from other devices may damage them.

Article configuration: Individually pouched.

In Chamber Conditioning:

Temperature:	55°C
Relative Humidity:	80%
Vacuum Set Point:	1.3 psia (9 kPa)

Preconditioning Time: 30 minutes

Exposure:

Temperature: 55°C
Relative Humidity: 80%
Ethylene Oxide Concentration: 740 mg/L
Gas Exposure Time: 60 minutes
Gas Type: 100% EO

Aeration:

Aeration Time: 24 hours
Aeration Temperature: 55 ± 4°C

STERRAD

Our recommendation with Sterrad is a ONE TIME sterilization. (If sterilized more than once the bio-material will become damaged and should not be used.) The device must be sterilized with a full cycle of the STERRAD 100S. The full cycle of the STERRAD 100S Sterilizer comprises TWO identical exposure phases (6 minute injection, 2 minute diffusion, and 2 minute plasma) preceded by a 10 minute pre-exposure plasma phase (1440 µL of 59% hydrogen peroxide). Follow the manufacturer's sterilization recommendations with the Sterrad 100S Sterilization System.

Implantation technique

If you are not familiar with this implant, you can use the technique suggested below, see how your experimented colleagues proceed, or use another type of

implant, because a faulty implantation technique can cause adverse effects. Although the choice of the implantation technique is the surgeon's responsibility, the recommendations below give important particular information about this implant.

In this procedure, both orbits are prepared and draped to allow for intraoperative comparison of volume. A subciliary or inferior forniceal approach, similar to that used for a blowout fracture, is developed. The incision is taken down to the periosteum, just anterior to the inferior orbital rim. The periosteum is incised and elevated, and the implant is slipped into the subperiosteal space (fig 2). Because of their configurations, a right implant will fit snugly behind the orbital rim in a right orbit. If a left implant is improperly selected, it will rock in place and fit poorly behind the rim.

With the implant in place, the ocular prosthesis is repositioned in the socket, the lid is replaced in its anatomic position, and both orbits are examined for volume symmetry. Generally, a slight overcorrection on the operating table is desired. If more volume augmentation is desired, the implant is removed and replaced with a larger one (larger bump). The implant is fixed with small titanium screws (7mm x 2.4mm) as it is simpler than wiring. The periosteum is closed, followed by skin or conjunctival closure, in the forniceal approach.

The addition of cold cure (like Cranioplast) acrylic in the O.R. is not advised by Oculo-Plastik, inc. That is why we offers three (3) sizes each side.

FIGURE 1 PREOPERATIVE APPEARANCE OF A PATIENT WITH A DEEP SUPERIOR SULCUS, FOLLOWING ENUCLEATION. POSTOPERATIVE APPEARANCE.


FIGURE 2 THE IMPLANT IS SLIPPED INTO THE SUBPERIOSTEAL SPACE.

PRECAUTIONS FOR USE

This implant is sold non sterile. Ethylene oxide sterilization (EO) is recommended.

The Codère-Durette implant is MR Safe and does not pose any known hazard in a magnetic resonance environment.

Codère-Durette is a **single use** medical device for a **single patient**:

- A second use of Codère-Durette at a later time may cause an infection due to microbial contamination.

The storage's duration is permanent.

Do not use any device that has been dropped accidentally without thorough examination of its condition. Do not use any device that is scratched or damaged.

Cold sterilization is not recommended for a bio-material intended for implantation.

The patient should visit yearly his ophthalmologist and his ocularist in order to minimize subsequent complications. If there are any changes in the performance of the implant, the patient must communicate with one or another of those specialists.

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ATTENTION, consulter les documents joints.

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