LASERSHIELD EYEWEAR

CAUTION, consult accompanying documents.

Get our instructions for use at oculoplastik.com.

oculoplastik.com

CAUTION:

Never look directly into the path of a laser beam. LaserShield eyewear offers protection against accidental exposure to stray or diffused reflection of laser beam energy for a maximum exposure of 10 seconds. LaserShield eyewear protects only the wearer. All personnel within the work area should wear appropriate eye protection against possible exposure to reflected beam energy. Increased work area lighting may be required if eye protection luminous transmittance is less than 20%.

All LaserShields are clearly marked with wavelength (nm) and absorption (OD) (Optical Density = $\log[1/T_{sig}]$ where T_{sig} =spectral transmittance) and L-ratings referencing tables 1 and 4 in EN207 (where applicable). Absorption curves available upon request. Color recognition (of warning lights, etc) may be impaired by use of tinted filters. LaserShield filters have passed environmental testing with extended temperature ranges of -65 to 160 degrees Fahrenheit.

CARE AND CLEANING:

- Store product in protective case when not in use.
- Store in areas not exceeding 80°F (26.6°C).
- Discard if damaged, faded or if scratches reduce vision.
- Clean with a mild detergent (e.g. Dawn) or any over the counter lens/sunglass cleaner (no alcohol) and wipe with non-abrasive cotton cloth.

PROTECTIONS VISUELLES LASERSHIELDS

ATTENTION, consulter les documents joints.

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LUMINOUS TRANSMITTANCE (%LT)
FOR EN207 COMPLIANT LASERSHIELDS

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