OCULOPLASTIK

(en) Durette implant 1

(fr) Implant Durette 14

Rxonly

(en) DURETTE IMPLANT

INSTRUCTIONS FOR USE

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

CAUTION, consult accompanying documents.


This Product is sold Non Sterile.

oculoplastik.com

Get our instructions for use at oculoplastik.com.

Do not autoclave.

Do not reuse.

Date of manufacture.

Rxonly

For prescription use only.

DESCRIPTION


The Durette is a PMMA ocular implant used for enucleation, evisceration, and secondary implantation (if small enough). It is available in two models, the Durette 1 and 4 in a 20 mm size only. (Fig A)

To distinguish how each anterior surface is designed:


- The Durette 1 is as spherical as possible. It has a slight double radius front (the horizontal being longer than the vertical) and is 1 mm flatter on the anterior surface compared to a sphere.
- The Durette 4 has still higher and larger mounds, although the mounds remain well within the sphere.

Its' basic design is quasi-spherical (Fig B). All models offer motility because of some frontal details for coupling. The permanently smooth surface of the implant should help decrease the possibility of potential erosion through covering tissue. The posterior has a medial off-center elongation to be positioned in the apex of the orbit or globe, to fill more and help placement of frontal details. This


elongation is slight so it does not interfere with motility.

They all have 3 sets of suturing tunnels. Situated in the anterior 3rd of the implant, they enable attaching the muscles to the ocular implant (Fig C). These 20 interconnected tunnels, between 1-2 mm in size, allow human tissue to grow inside, help with stability and limit the risks of migration and extrusion. Note that we had to make intermediate level sets different in model 4 compared to model 1.

The size of any implant should be chosen, given available space in the orbit, and as long as closing of the tissues is possible without tension. Although the surgeon makes the choice of model and size, we now have only the 20 mm available in the 2 models. If impossible to close tissue without tension, another smaller implant must be used. The front details should allow sufficient space for the prosthesis with as normal an anterior chamber as possible.

APPEARANCE

The Durette implant is delivered non sterile. It should not be cleaned prior to sterilization.

ACCESSORIES

- Synthetic absorbable Polysyn™ braided coated polyglycolic acid, violet, 18"/45 cm, 5-0 double-armed needle HRT17 taper-cutting ½ circle 17 mm (Sharp point #GA-502N (Angiotech). – not included. These essential needles pass well in all tunnels. Angiotech accepted to make this special needle/suture combination for this implant.
- Four Herrera clips to hold the implant during surgery, 12-020 – not included.
- Acrylic Conformer 22 (small or regular) for all models, 19-092 or 19-094 – not included.

INDICATIONS

The Durette ocular acrylic (PMMA) implants in 2 models are permanent implants that occupy the eye cavity when it becomes necessary to surgically remove the eye (enucleation), the contents of the eye sac (evisceration), or space left after the removal of another ocular implant (used as a secondary implant). It is used to replace volume and to impart motion and stability to the eventual ocular prosthesis.

MODE OF ACTION

Both models offer increased motility because of some frontal details for coupling. The posterior medial off-center elongation must be positioned in the apex, to help placement of frontal details, fill more space and push the orbital intraconal fat posteriorly and superiorly. They all have suturing tunnels to attach the muscles to the implant. The suggested HRT17 needle works best for all the tunnels. The network of interconnected tunnels allows human tissue to grow inside and help with stability. Once incisions in posterior Tenon's are done, there is sufficient tissue needed to cover the implant for a closure without tension.

CONTRAINDICATIONS

Contraindications for this implant are: cases of severe scarring (pemphigoid, trachoma, and injuries), in which the patient has insufficient Tenon's capsule and/or conjunctiva to properly cover the implant. Only smaller implants should be used in such cases, as sufficient tissue is needed to cover the implant for a closure without tension.

In parts of the world where ocular prostheses are not usually custom-made, only the Durette-1 implant is advised. In such cases, ideally the back of the ocular prosthesis must be adjusted to avoid

any pressure points. Fitting as if it were a sphere could be done although it is still recommended to custom impression-fit the ocular prosthesis.

ADVERSE EFFECTS

As with other PMMA implants, complications should be limited and the attachment of the muscles to the implant should reduce migration, which is often a problem with smooth surface conventional spherical implants. Complications can still include and not be limited to:

- Tissue dehiscence or erosion
- Implant exposure
- Implant extrusion
- Deep superior sulcus deformity
- Volume deficiencies
- Socket cul-de-sacs contraction
- Ptosis of the superior and inferior eyelid
- Poor motility
- Migration or displacement of the implant
- Discomfort
- Infection

The risk of extrusion for this implant is expected to be very low, because of its permanent smooth surface anteriorly.

New tissue invasion in the interconnected tunnels may also decrease the extrusion rate. These ingrown tissues stabilize it in place, minimizing migration, including forward migration and stretching of the covering tissues. Undue tension of the covering tissue plays a role in extrusion.

It is recommended that a custom impression-fitted ocular prosthesis should be made over both models of the Durette implant. If a custom fit prosthesis cannot be made, a Durette 1 should be used, unless the patient is able to travel to obtain a custom-made impression-fitted ocular prosthesis. Unexpected side effects and complications related to this medical device must be reported to Oculo-Plastik, Inc. If you have any suggestions to improve these Instructions For Use, please let us know.

RECOMMENDATIONS FOR USE

STERILIZATION

Follow your specific sterilizer equipment manufacturer's recommendations, or use cycles that have been previously established and validated for your facility.

The Non Sterile device should remain in the original sealed packaging until it will be STERILIZED. It should not be cleaned prior to sterilization. It

should be sterilized alone in a separate Peel Pouch that does not contain other devices.

- The device can only be sterilized one time.
- Any device that has become contaminated by blood, tissue and/or bodily fluids should be thrown out and should not be processed in any way.
- The device should not be processed or transported or come in contact with any type of contaminated devices and/or equipment.
- The device should be transported separate from soiled devices in its original sealed packaging in order to avoid contamination.
- The device should remain in the original packaging, which is sealed until it will be STERILIZED in order to avoid contamination.
- The sealed package should be opened under clean conditions.
- Care must be taken to prevent possible contamination by talc, dust, lint, surgical glove powder, drape and sponge lint, fingerprints, skin oils, and other surface contaminants that may be deposited during handling prior to sterilization.
- Talc, dust, lint, surgical glove powder, drape and sponge lint, fingerprints, skin oils and other surface contaminants deposited by improper

handling may adversely affect the suitability of the device.

- Strict adherence to clean, aseptic techniques should be maintained to prevent contamination of the device and possible complications.
- ONCE THE DEVICE HAS BEEN REMOVED FROM THE ORIGINAL PACKAGING IT SHOULD BE STERILIZED IMMEDIATELY in a SEPARATE peel pouch that does not contain other devices.

ETHYLENE OXIDE

To sterilize pouched devices, place them on the sterilizer rack, in a single layer, as pressure from other devices may damage them.

Article configuration: Individually pouched.

In Chamber Conditioning:

Temperature: 55°C
Relative Humidity: 80%
Vacuum Set Point: 1.3 psia (9 kPa)
Preconditioning Time: 30 minutes

Exposure:

Temperature: 55°C
Relative Humidity: 80%
Ethylene Oxide Concentration: 740 mg/L

Gas Exposure Time: 60 minutes
Gas Type: 100% EO

Aeration:

Aeration Time: 24 hours
Aeration Temperature: $55 \pm 4^{\circ}\text{C}$

Implantation technique

The technique suggested below is recommended and gives information specific to this implant.

ENUCLEATION

A 360° peritomy is done in a sub Tenon plane down to the sclera as close as possible to the cornea. The Stevens scissors are used to gently spread the tissue in each quadrant and to sever any remaining adhesions between the globe and Tenon’s capsule. After Tenon’s capsule has been freed from the muscle insertions over 5 mm, the recti muscles are mounted on 5-0 Polysyn™ (polyglycolic acid) suture from Sharpoint Angiotech (45 cm with double HRT17 needles, item #GA-502N). They are then cut at their insertions. No further dissection of Tenon’s capsule from the muscle should be done. All rectus muscles are severed in this manner. It is very im-

portant that Tenon’s capsule and the intermuscular septum are not cleaned from the muscle bellies.

The inferior oblique muscle and superior oblique tendon are detached from the globe before its removal. A 1 to 2 mm stump of the medial rectus muscle tendon may be left on the globe to facilitate grasping the globe for cutting the optic nerve stump. A suture passed into the insertion site of one of the medial rectus muscles may also help.

Removal of the Globe

The globe is gently pulled anteriorly. The scissors are then placed posteriorly to the globe and used to spread the tissue on either side of the optic nerve. This maneuver frees the posterior pole of the globe. The closed scissors are used to determine the position of the optic nerve. The scissors are then opened, one blade placed on either side of the optic nerve, and the optic nerve is severed. Some surgeons may prefer clamping the nerve before cutting it. Any remaining tissue attachments are severed as close to the globe as possible as it is delivered.

Tamponade is achieved with gauze held firmly in the enucleation cavity for approximately five minutes, or until bleeding stops. Make sure that

the posterior off-center elongation is in the apex of the orbit. To allow placement of the 20 mm implant and allow ample closure of anterior Tenon's, posterior Tenon's is opened, as necessary, to minimize tension of the tissues covering the implant. A 20 mm sizer may be used to determine if this 20 mm Durette implant will be used.

Precautions

The success of the implant is increased by the following precautions in removing the globe:

1. Avoid dissection of Tenon's capsule from the muscles and keep the intermuscular septum intact;
2. All active bleeding must be stopped by pressure application with a tamponade wrapped in gauze, since bleeding within Tenon's capsule after implantation increases the possibility of extrusion; and
3. Incisions in posterior Tenon's must be made to obtain proper closure over the chosen implant.

Insert a 20 mm stainless steel sizer (Oculo-Plastik) in the orbit to validate proper choice of size of the 20 mm Durette implant.

Optional suturing platform

During surgery, surgeons may elect to first making a suturing platform on the implant while working on the sterile instruments table in the O.R. The HRT17 needle enables easily passing the Polysyn™ absorbable sutures in the chosen sets of holes of the implant (Figure F). Care must be taken that the suture platform has multiple separately knotted sutures. Individually knotted sutures are pulled so the knot is hidden inside the tunnel. The platform can be in a different configuration as preferred. The platform helps avoid having to suture the muscles directly to the implant tunnels. Instead surgeons will suture muscles to sutures already on the platform. Making a platform is an option. For safety reasons, we did not make it in advance. Surgeons do not have to make a suturing platform as the HRT17 needle passes well inside the tunnels. Simply follow the ceiling of the tunnels to facilitate passage.


FIGURE F

Durette ocular implant - Before and after making suture platform (knots of individual sutures are hidden inside tunnels).

The implant should be soaked in an antibiotic solution prior to implantation.

Proper medial positioning of posterior elongation with the help of the Herrera holding clips

Whether they make a suturing platform or not, surgeons should first establish proper implant orientation, with the posterior off-center elongation being medial. One Herrera clip should be installed to maintain proper implant orientation during surgery and to insure stability while suturing (Figure G).

FIGURE G

Durette ocular implant - With one Herrera clip to hold implant while validating closure.

FIGURE H

Durette ocular implant - photographed with the four Herrera clips in each quadrants to hold the implant steady during surgery and maintain proper orientation.

Once with one clip in place, the implant is inserted deep inside the socket, to again validate ample closure. More than one Herrera clip (Figure H) can then be inserted in the lateral holes in the quadrants. The clips can be removed easily.

Muscles positioning if a platform is made

If a suturing platform is made, muscle positioning will depend on the surgeons' platform and configuration.

One arm of each double-arm 5-0 Polysyn™ suture with HRT17 needle (Angiotech Sharpoint, item #GA-502N) is secured to the suturing platform on the implant for each of the four rectus muscle and each muscle is secured to each other. More details follow.

Muscles positioning if sutured directly to implant tunnels

If muscles are attached directly to the implant, one arm of each double-arm 5-0 Polysyn™ suture with HRT17 needle is passed in chosen tunnels of the implant and each rectus muscle is secured to the implant and to each other (as shown in the figure I).

All the muscles are attached in the central position with a square knot. If the implant is seen to tilt to

one side or the other, the muscles can be reattached to compensate before the permanent knots are completed.

Once this is done, if desired, to further secure the muscles to the implant and to each other, additional sutures can be passed in the holes of the implant.

The lateral tunnels can be used for cases where a muscle has been shortened or scarred (from trauma, incarceration in a fracture site or from previous surgery).



Figure I shows the suggested muscles layout with sutures from tunnels placed through from belly surface of rectus muscles. Vertical muscles are attached first, then horizontal. Note that the verti-

cals are sutured to the central set of tunnels without touching while the horizontals, also sutured to the central set of tunnels, are touching.

Some surgeons will prefer tying the superior rectus to the lateral set of tunnels to avoid pulling the levator and avoid causing or worsening a possible ptosis of the upper lid. (Figure J)


Others, who feel it may offer better support, may also put a double armed suture through the inferior oblique muscle stump after cutting it off and tack it to the inferior edge of the lower rectus. (Not illustrated)

With additional interrupted sutures from close-by tunnels, the horizontal to the vertical rectus muscle can be sutured to each other in each quadrant,

if stability of the implant is warranted further (Figure K).

Suggested final muscles layout with spread out muscles in each quadrant and additional suturing of muscles to intermediate holes. Note the inevitable difference in suturing the intermediate tunnels for the model 1 versus the model 4.

Closure

The proper and complete closure of Tenon's capsule is the most important factor in the permanence of any implant. Tenon's capsule can be closed quickly and completely with multiple, interrupted 4-0 or 5-0 Synthetic absorbable Polysyn™ (polyglycolic acid) suture in the horizontal plane. Meticulous closure of Tenon's is important and

should not allow any visualization of the implant. Conjunctiva is closed separately. Opening in Tenon's capsule are easily found by direct observation or by probing the closure sites with the tip of a muscle hook. The conjunctiva is closed horizontally with a running 6-0 plain suture. Some prefer meticulously closing both Tenon and conjunctiva with multiple interrupted sutures. Closure other than horizontal reduces the width of the cul-de-sac and limits the horizontal movement of the prosthesis. The size of the implant chosen should allow closure of Tenon's and conjunctiva without tension.

An antibiotic ointment is instilled followed by the placement of a conformer. An intermarginal suture is recommended in case a postoperative hemorrhage happens. The patient is released after adequate pain control has been achieved and is seen 2 to 3 days after surgery to rule out infection and remove the intermarginal suture and then at 1 week and 6 weeks postoperatively. A custom fit prosthesis is fabricated when healing is adequate (at 12 weeks). A custom fit prosthesis is fabricated when healing is adequate.

EVIScerATION

A 360° limbal peritomy is made using Westcott scissors and staying as close as possible to the limbus. Then, the conjunctiva is dissected 5-10 mm away from limbus. A number 11 blade is inserted in the limbus to enter the anterior chamber, and a 360° keratectomy is carried out with Westcott's scissors. An evisceration spoon is inserted in the suprachoroidal space and using semicircular movements is directed more posteriorly around uvea and vitreous body. An attempt is made to remove the contents of the eye in one piece although this is not always possible. Once the major portion of the intrascleral contents has been removed, the bleeding vessels are visualized and cauterized if needed. Care is taken to remove all residual uveal tissue.

In cases of endophthalmitis, active infection should be treated before implantation is performed. Cautery is insured throughout the surgery.

Posterior sclerotomies are created all around the optic nerve to disinsert it. Radial sclerotomies are also created around the optic nerve stump. Two anterior sclerotomies are made at 3 and 9 o'clock at the limbus to facilitate implant placement. Spherical sizers may help in choosing the right size of the definitive implant. Once the sclera can be

overlapped without tension, the implant is secured with Polysyn™ 5-0 (polyglycolic acid) (45 cm, single arm, with HRT17 needles to be able to penetrate tunnels). A suture platform can also be made during evisceration surgery to help secure the implant with details in front. While fixing the Durette implant in the desired position, insure medial off-center elongation is in the apex and the front details well in the front. A partial thickness scleral bite is taken superiorly and inferiorly ~8 mm inside the scleral lip with the pre-placed sutures sitting in the superior and inferior intermediate set of tunnels. When tied, these sutures will hold the implant in position until sclera contracts down on it. Relaxing incisions with removal of “dog ears” is often necessary at 3:00 and 9:00 o’clock. Three mattress sutures with Polysyn™ 5-0 (polyglycolic acid) are used to close sclera anteriorly with a few millimeters of overlap. Closure is completed with multiple synthetic absorbable Polysyn™ 4-0 as simple stitches to reinforce the overlapped scleral closure. When the sclera is closed in this fashion, there will still be some space between the Durette implant and sclera, but, over the ensuing 4-8 weeks, the sclera will contract and the space between sclera and implant will disappear. Within 2 months, the sclera shell will be wrapped tightly against the Durette implant.

The conjunctiva is now closed with a 6-0 plain suture being sure to cover the synthetic absorbable Polysyn™ (polyglycolic acid) suture ends with conjunctiva. An antibiotic ointment is instilled followed by the conformer and an intermarginal suture if needed. More topical antibiotic ointment is instilled in the eye and a compressive patch is applied. The patient is released after adequate pain control has been achieved and is seen 2 to 3 days after surgery to rule out infection and remove the intermarginal suture if present and then at 1 week and 6 weeks postoperatively. A custom fit prosthesis is fabricated when healing is adequate.

An acrylic conformer should be used with this implant. The use of a silicone conformer is never recommended.

PRECAUTIONS FOR USE

This implant is sold non sterile. Ethylene oxide sterilization (EO) is recommended. Do not autoclave.

The Durette implant is MR Safe and does not pose any known hazard in a magnetic resonance environment.

The Durette implant is a **single use** medical device for a **single patient**:

- A second use of Durette implant at a later time may cause an infection due to microbial contamination.

The storage's duration is permanent.

Do not use any device that is damaged.

Cold sterilization is not recommended for a bio-material intended for implantation.

The patient should visit yearly his ophthalmologist and his ocularist in order to minimize subsequent complications. If there are any changes in the performance of the implant, the patient must communicate with one or another of those specialists.

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ATTENTION, consulter les documents joints.

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