OCULOPLASTIK

oculoplastik.com

(en) THERMOFORMING SHEET

INSTRUCTIONS FOR USE

READ ALL INSTRUCTIONS BEFORE USING THIS PRODUCT.

CAUTION, consult accompanying documents.

oculoplastik.com

Get our instructions for use at oculoplastik.com.

Do not reuse.

Description

Cell cast acrylic is a rigid thermoplastic material which is characterized by its lit-

tle aging and discoloration and has excellent UV stability. It is used by the ocularists to make a fitting shell (transparent sheet) or could be part of a ultra-thin ocular prosthesis (white sheet).

Appearance

Available in transparent or white 2 x 12 x 1/16 inches (50 x 305 x 1.5 mm) sheets.

Indications

Depending on the sensitivity of the eye being fitted, the ocularist must decide if a clear trial shell must be made or if the prosthesis must be so thin that the only method possible to make the prosthesis is the one using the opaque white acrylic material to press it. The clear shell has the major advantage of providing a window through which the quality of fit and the reaction of the globe can be visualized. Going directly from the impression to a white shell eliminates an additional duplication of the shell, which could result in the loss of detail.

Mode of action

Heated to a pliable state, pressed, cooled, trimmed and polished. When hot and flexible the thermoforming sheet is pressed onto the dental stone mould of the patient's globe with a scleral shell press. The obtained shape is trimmed, polished and washed before being worn by the patient.

Contraindications

Not for use over seeing healthy eyes. This material is to make devices over shrunken or not, unsightly or not, non-seeing globes or a globe with no useful vision.

Advice for patient

These forms must never be left alone under high heat or the sun, as they may distort and regain their flat or near flat original state.

Adverse effects

Proper fit depends on the work of the ocularist. Improperly made shell or

prosthesis may press, hurt, damage the globe of the patient.

Recommendations for use

Once the impression of the globe is taken, dental stone is poured into a small cup and the impression is cast to make the model of the globe. An electric oven range is used to heat the end of the thermoforming sheet and once the sheet is pliable, rapidly, the soft sheet is placed over the stone model of the globe and the rubber of the press is lowered and locked so that the sheet covers the stone completely. Once cooled, the moulded sheet is removed from the stone and verified that it did cover the stone replica totally. The form is trimmed and polished. Both the clear and white sheets are moulded the same way. The clear is to become a trial shell over a shrunken or not globe, non-seeing globe or a globe with no useful vision. The white is to become the shape of the ocular prosthesis. The white is thinned in front. The iris, sclera tinting and veining is done or painted on and as usual in ophthalmology, clear polymer

acrylic dough is cured over the artwork to seal it.

Precautions for use

Polish well. Wash well before inserting. Be careful when placing in patient eye socket.

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ATTENTION, consulter les documents joints.

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OCULO-PLASTIK, Inc.